

Lifeloc WorkPlace Series



Lifeloc

Professional Breath Alcohol Tester

With the addition of Bluetooth capability for the Phoenix® 6.0BT evidential breath testing system, you now have access to secure wireless printing with unparalleled mobility and ease of use.

A fully integrated solution, the Phoenix 6.0BT comes standard with Lifeloc's proprietary EasyMode™ software, guiding users through the DOT testing protocol, and our PermAffix Labeling System.™ The 6.0BT is the most advanced and powerful evidential breath tester on the market.

Phoenix 6.0BT



Unlock the Power of Alcohol Testing

Lifeloc WorkPlace Series



PHOENIX 6.0BT

Lifeloc has paired our signature technologies, EasyMode Software and PermAffix Labeling with Bluetooth wireless printing, providing the ultimate DOT and Workplace testing solution.

EasyMode Software

- **Ease of Use:**
Provides error free testing by guiding you through DOT protocol. (Automatic Air Blank - Screening Test - 15 Minute Timer - Confirmation Test - Printout)
- **Ensures Compliance:**
On a positive test result, prompts are given when you need to contact the DER as well as when to perform a calibration check

PermAffix Labeling System

- **Tamper Evident:**
Slits on all four corners providing enhanced security reducing the need for messy tape
- **Self Adhesive:**
Simply peel and stick with a powerful adhesive protecting against peeling
- **Sized to Fit:**
Cut specifically to fit DOT and other breath alcohol testing forms
- **Permanent:**
Clear imaging without smudging or fading

Wireless Bluetooth Printing

- **True Mobility:**
Up to 33' wireless connectivity range
- **Secure:**
Encrypted data transmission between printer and unit
- **Automatic Temperature Compensation:**
Heated print head provides fast, clear results in any climate
- **Long Battery Life:**
Low power draw from the printer
- **Automatic Bluetooth Pairing:**
Your 6.0BT is paired with your printer before leaving the factory. Simply turn the unit on and you are ready to print.

Other Features

- **Deluxe Warranty**
- **Automatic Altitude Compensation**
Built in barometric pressure sensor for easy calibration
- **Bright Graphic Display**
- **Fast Response and Recovery Times**
- **Customizable Software**
- **Platinum Fuel Cell Sensor**
- **Five Star Customer Care**
Unbeatable service, all repairs are 5 days or less or service is free
- **Training**
National network of Certified Master Trainers to meet your specific needs
- **DOT/NHTSA Approved**
- **AlcoMark® Software Compatible**

